



Saybolt Colorimeter LD-LSBO-A10

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Saybolt Colorimeter LD-LSBO-A10

Overview

Saybolt Colorimeter LD-LSBO-A10 offers a measurement range from -16 (deepest) to +30 (lightest). Features color numbers with corresponding liquid column heights for accurate color assessment. Includes both half-thick and full-thick standard color plates on a rotating table.

Features:

- Equipped with a graduated tube
- Securely Mounted Plain Tube
- Includes both half and full-thickness standard swatches
- Exquisite structure and excellent production
- Operation is adaptable and convenient
- Test results are quick, clear, and accurate

Applications:

Saybolt Colorimeter is widely used in petroleum and chemical industries essential for quality control of petroleum products like crude oil, gasoline, diesel, and lubricants.

Specifications:

Weight	12 Kg
Dimensions	140 × 140 × 670 mm
Power Supply	AC 220-240 V, 50/60 Hz
Working Mode	Standard Color Tube and Sample Glass Tube with Scale
Light Source Lamp	5V/2A
Measurement Range	-16 (deepest) to +30 (lightest)
Power Consumption	100 W
Compliant with Standards	GB/T3555, ASTM D156