



## **Polypropylene Dropper Bottle LD-PDB101L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polypropylene Dropper Bottle LD-PDB101L

## Features:

- Material: Polypropylene and Low Density Polyethylene
- Leak proof sealing design
- Pyrogenic free
- Autoclavable

## Specifications:

<b>Drop</b>	50 ?L
<b>Height</b>	60
<b>Volume</b>	8 mL
<b>Diameter</b>	12.5
<b>Packaging</b>	Plastic Bag
<b>Bottle Color</b>	Transparent / White
<b>Quantity per Case</b>	3000