



Polypropylene Dropper Bottle LD-PDB104L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Polypropylene Dropper Bottle LD-PDB104L

Features:

- Material: Polypropylene and Low Density Polyethylene
- Leak proof sealing design
- Pyrogenic free
- Autoclavable

Specifications:

| | |
|--------------------------|---------------------|
| Drop | 50 ?L |
| Height | 61 |
| Volume | 8 mL |
| Diameter | 20 |
| Packaging | Plastic Bag |
| Bottle Color | Transparent / White |
| Quantity per Case | 4000 |