



**Polyethylene Square PET Bottle LD-
PSB101L**

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polyethylene Square PET Bottle LD-PSB101L

Features:

- Material: USP Class VI Polyethylene Terephthalate
- Screw Cap Material: High Density polyethylene
- Sterilization Level is 10^{-6}
- Light weight and shatter proof
- Clear graduation for accurate measurement
- It has low permeability to carbon dioxide and oxygen
- DNase, RNase and pyrogenic free

Specifications:

Height	147 mm
Capacity	250 mL
Packaging	32 pcs/Bag
Quantity per Case	96