



**Polyethylene Square PET Bottle LD-  
PSB103L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polyethylene Square PET Bottle LD-PSB103L

## Features:

- Material: USP Class VI Polyethylene Terephthalate
- Screw Cap Material: High Density polyethylene
- Sterilization Level is  $10^{-6}$
- Light weight and shatter proof
- Clear graduation for accurate measurement
- It has low permeability to carbon dioxide and oxygen
- DNase, RNase and pyrogenic free

## Specifications:

<b>Height</b>	208 mm
<b>Capacity</b>	1000 mL
<b>Packaging</b>	24 pcs/Bag
<b>Quantity per Case</b>	48