



Polyethylene Wash Bottles LD-PWB100L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polyethylene Wash Bottles LD-PWB100L

Features:

- Material: Low Density Polyethylene
- Tip can be cut back to provide increased flow
- Colored Screw Cap (Green, Blue, Red, Yellow, White)
- Chemically resistant to strong acids or alkali

Specifications:

Height	110 mm
Capacity	250 mL
Mouth Diameter	0.8 mm
Quantity Per Case	50