



## **Polyethylene Wash Bottles LD-PWB102L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polyethylene Wash Bottles LD-PWB102L

## Features:

- Material: Low Density Polyethylene
- Tip can be cut back to provide increased flow
- Colored Screw Cap (Green, Blue, Red, Yellow, White)
- Chemically resistant to strong acids or alkali

## Specifications:

<b>Height</b>	170 mm
<b>Capacity</b>	1000 mL
<b>Mouth Diameter</b>	0.8 mm
<b>Quantity Per Case</b>	50