



Plastic Sample Container LD-PSC205L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Plastic Sample Container LD-PSC205L

Features:

- Material: Polystyrene or Polypropylene
- Cap Material: Polypropylene
- Exclusively for In Vitro Uses
- Leak-proof Tested
- EO Sterilized

Specifications:

Type	Plain Label
Capacity	30 mL
Material	Polypropylene
Packaging	50 Pcs/Bag, 8 Bags/Case
Quantity Per Case	400