



## **Plastic Sample Container LD-PSC305L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Plastic Sample Container LD-PSC305L

## Features:

- Material: Polystyrene or Polypropylene
- Cap Material: Polypropylene
- Exclusively for In Vitro Uses
- Leak-proof Tested
- EO Sterilized

## Specifications:

<b>Type</b>	Plain Label
<b>Capacity</b>	60 mL
<b>Material</b>	Polypropylene
<b>Packaging</b>	60 pcs/Tray, 5 Tray/Case
<b>Quantity Per Case</b>	300