



## **Plastic Test Tube LD-PTT100L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Plastic Test Tube LD-PTT100L

## Features:

- Chemically resistant
- Attainable maximum RCF 5000 x g without any crack or leakage
- Suitable for storage of organic solvents, weak acids and bases
- Labels can be customized

## Specifications:

<b>Capacity</b>	5.5 mL
<b>Material</b>	Polypropylene
<b>Packaging</b>	400 pcs/bag, 10 bags/case
<b>Quantity Per Case</b>	4000
<b>Tube Dimension ( ? x h )</b>	13 x 78 mm