



Plastic Test Tube LD-PTT103L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Plastic Test Tube LD-PTT103L

Features:

- Chemically resistant
- Attainable maximum RCF 5000 x g without any crack or leakage
- Suitable for storage of organic solvents, weak acids and bases
- Labels can be customized

Specifications:

Capacity	8 mL
Material	Polypropylene
Packaging	250 pcs/bag, 12 bags/case
Quantity Per Case	3000
Tube Dimension (? x h)	15 x 100 mm