



Plastic Test Tube LD-PTT110L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Plastic Test Tube LD-PTT110L

Features:

- Chemically resistant
- Attainable maximum RCF 5000 x g without any crack or leakage
- Suitable for storage of organic solvents, weak acids and bases
- Labels can be customized

Specifications:

Capacity	10 mL
Material	Polystyrene
Packaging	250 pcs/bag, 8 bags/case
Quantity Per Case	2000
Tube Dimension (? x h)	15 x 100 mm