



Polypropylene Beaker LD-PBK102L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polypropylene Beaker LD-PBK102L

Features:

- Material: Polypropylene
- Autoclavable at 120 °C
- Rigid and translucent, enabling easy and consistent viewing of liquid levels
- Chemically Resistant
- Graduated in mL

Specifications:

Capacity	50 mL
Packaging	100 pcs/Bag, 10 Bags/Case
Graduation Range	5 - 50 mL
Quantity Per Case	1000