

# **Polypropylene Beaker LD-PBK104L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polypropylene Beaker LD-PBK104L

## Features:

- Material: Polypropylene
- Autoclavable at 120 °C
- Rigid and translucent, enabling easy and consistent viewing of liquid levels
- Chemically Resistant
- Graduated in mL

## Specifications:

<b>Capacity</b>	150 mL
<b>Packaging</b>	100 pcs/Bag, 10 Bags/Case
<b>Graduation Range</b>	25 - 150 mL
<b>Quantity Per Case</b>	1000