



Polypropylene Beaker LD-PBK105L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Polypropylene Beaker LD-PBK105L

Features:

- Material: Polypropylene
- Autoclavable at 120 °C
- Rigid and translucent, enabling easy and consistent viewing of liquid levels
- Chemically Resistant
- Graduated in mL

Specifications:

Capacity	200 mL
Packaging	40 pcs/Bag, 10 Bags/Case
Graduation Range	50 - 200 mL
Quantity Per Case	400