



Polypropylene Beaker LD-PBK106L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Polypropylene Beaker LD-PBK106L

Features:

- Material: Polypropylene
- Autoclavable at 120 °C
- Rigid and translucent, enabling easy and consistent viewing of liquid levels
- Chemically Resistant
- Graduated in mL

Specifications:

Capacity	250 mL
Packaging	50 pcs/Bag, 10 Bags/Case
Graduation Range	25 - 150 mL
Quantity Per Case	500