



Polystyrene Petri Dish LD-PPD103L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polystyrene Petri Dish LD-PPD103L

Features:

- Material: Polystyrene
- Serrated ring for easy handling
- Number codes for easy identification
- Inner surface molded with vent for air exchange
- Pre-sterilized
- DNase, RNase, Pyogenic and toxic free

Specifications:

Dimension	100 x 20 mm
Packaging	300 pcs/bag
Culture Area	55.8 cm ²
Quantity per case	3000