



## **Polystyrene Petri Dish LD-PPD104L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polystyrene Petri Dish LD-PPD104L

## Features:

- Material: Polystyrene
- Serrated ring for easy handling
- Number codes for easy identification
- Inner surface molded with vent for air exchange
- Pre-sterilized
- DNase, RNase, Pyogenic and toxic free

## Specifications:

<b>Dimension</b>	90 x 16 mm
<b>Packaging</b>	500 pcs/bag
<b>Culture Area</b>	53.2 cm <sup>2</sup>
<b>Quantity per case</b>	5000