



Polystyrene Solution Reservoirs LD-PSR106L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polystyrene Solution Reservoirs LD-PSR106L

Features:

- Material: Virgin Polystyrene
- Designed for use with 8-channel and 12-channel pipette
- Equipped with poring spouts in each corner
- Tapered bottom providing a narrow tapered trough to minimize reagent loss and maximize recovery
- Available in Non-sterilized and Electron Beam Sterilized

Specifications:

Capacity	50 mL
Packaging	100 pcs/Case
Sterility	Non-Sterile
Quantity Per Case	100