



Polystyrene Weighing Boats LD-PSW102L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polystyrene Weighing Boats LD-PSW102L

Features:

- Material: Polystyrene
- Ideal for handling both Solids and liquids
- Rounded edges for easy removal of samples
- Available in white and black color

Specifications:

Type	Square Shaped
Color	White
Volume	250 mL
Dimension	140 x 140 mm
Packaging	250 Pcs/Bag, 4 Bags/Box, 4 Box/Case
Quantity Per Case	4000