



## **Sampling Bag LD-SBG111L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Sampling Bag LD-SBG111L

## Features:

- Material: Polyethylene
- Resistant to high temperature
- Sterilized by gamma irradiation
- DNase, RNase and Pyrogen free

## Specifications:

<b>Size</b>	300 x 180 mm
<b>Capacity</b>	1650 mL
<b>Packaging</b>	250 pcs/Box, 4 Box/Case
<b>Thickness</b>	106 ?m
<b>Writing area</b>	No
<b>Wire Thickness</b>	0.5 mm
<b>Quantity per Case</b>	1000