



Sampling Bag LD-SBG200L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Sampling Bag LD-SBG200L

Features:

- Material: Polyethylene
- Resistant to high temperature
- Sterilized by gamma irradiation
- DNase, RNase and Pyrogen free

Specifications:

Size	127 x 76 mm
Capacity	60 mL
Packaging	500 pcs/Box, 2 Box/Case
Thickness	70 ?m
Writing area	Yes
Wire Thickness	0.5 mm
Quantity per Case	1000