



Sampling Bag LD-SBG201L

Labodam Equipment Ltd.

info@labodam.com | www.labodam.com

Sampling Bag LD-SBG201L

Features:

- Material: Polyethylene
- Resistant to high temperature
- Sterilized by gamma irradiation
- DNase, RNase and Pyrogen free

Specifications:

Size	178 × 76 mm
Capacity	120 mL
Packaging	500 pcs/Box, 2 Box/Case
Thickness	70 ?m
Writing area	Yes
Wire Thickness	0.5 mm
Quantity per Case	1000