



## **Vertical Laboratory Autoclave LD-LVA-K10**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Vertical Laboratory Autoclave LD-LVA-K10

## Overview

Vertical Autoclave LVA-K10 is a top loading sterilizer which provides safe, economic and effective sterilization of laboratory instruments. The autoclave has automatic system for feeding water, pulse vacuum, sterilization, vacuum drying, exhaustion and draining. These autoclave comes with digital display with touch key pad.

## Features:

- 7 inches' touch screen display
- Automatic control system
- Elevated type door
- Automatic system for feeding water, pulse vacuum, sterilization, vacuum drying, exhaustion and draining
- Equipped with PT-TT testing port
- Energy efficient steam generator
- Equipped with 0.2 $\mu$ m sterile ventilation filter
- Clean and dry sterilizing environment facilitated by steam-water inner circulation system

## Applications:

This product finds applications in: Biochemistry, Pharmaceutical Industry, Medicine and Health care Industry, Biotechnology, Microbiology and Scientific Research.

## Specifications:

<b>Capacity</b>	60 L
<b>Timer Range</b>	0 - 99 min
<b>Heat Average</b>	? $\pm$ 1 ?
<b>Sterilization Pressure</b>	0.22 MPa
<b>Sterilization Temperature</b>	134 °C
<b>Chamber Dimension ( ? x D )</b>	396 x 490 mm