



## **Rapid Low Temperature Circulating Bath LD-LRTB-A40**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Rapid Low Temperature Circulating Bath LD-LRTB-A40

## Overview

Rapid Low Temperature Circulating Bath LRTB-A40 is a microcomputer controlled 7 L bath with refrigeration system and a temperature range of - 40 °C to 200 °C. Equipped with dual red and green LED windows with upper window measuring red value and lower measuring green value. Features internal and external cycles with external cycles discharging thermostatic liquid outside the bath raising the temperature of the external vessel.

## Features:

- Volume - 7 L
- Low temperature cooling pump
- Microcomputer controlled temperature
- Temperature and time LCD display
- Built-in liquid level monitor
- Refrigeration unit with relay , protector and capacitor
- Anti-corrosive circulatory system

## Specifications:

<b>Volume</b>	7 L
<b>Pump flow</b>	20 L / min
<b>Heating power</b>	1200 W
<b>Head of delivery</b>	10 m
<b>Temperature range</b>	- 40 °C to 200 °C
<b>Temperature stability</b>	± 0.1 °C
<b>Refrigerating capacity</b>	2329 W
<b>Temperature control tier number</b>	32