



## **Walk-in Fumehood LD-LWF-A11**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Walk-in Fumehood LD-LWF-A11

## Overview

Walk-In Fumehood LWF-A11 has a large opening and is mounted on the floor. It is ideal in situations that require large materials to be walked or loaded in. It can be opened to a maximum height of 1300 mm. Meets the standard of a low velocity high performance hood. Consists of large room exhaust system and option of filter installation which minimizes environmental pollution.

## Features:

- 2 unit of LED lamp: 12 W ~ 16 W
- Memory function in case of power failure
- Operators can walk into the fumehood
- Large room exhaust system
- Option of filter installation which minimizes environmental pollution
- Manual toughened front glass window
- Exterior material made up of 1.0 mm cold-rolled steel
- Interior made up of excellent quality melamine board resistant to acid and base

## Applications:

Walk-in Fumehood find applications in pharmaceutical, educational, hospital, healthcare and manufacturing laboratories.

## Specifications:

<b>Noise</b>	62 dB
<b>LED Lamp</b>	16 W × 2
<b>Vibration</b>	XYZ direct ? 4 µm
<b>Illumination</b>	? 800 Lux
<b>Maximum Opening</b>	1300 mm
<b>Power Consumption</b>	600 W
<b>Exhaust Blower Size</b>	970 × 810 × 720 mm
<b>Inflow Air Velocity</b>	0.3 m/s to 0.8 m/s
<b>Exhaust Airflow Volume</b>	1500 m <sup>3</sup> /h
<b>Package Size (W × D × H)</b>	1710 × 1050 × 2430 mm
<b>Exterior Size (W × D × H)</b>	1500 × 800 × 2200 mm
<b>Interior Size (W × D × H)</b>	1280 × 600 × 1650 mm