



## **Mini Handheld Homogenizer LD-LMHH-A10**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Mini Handheld Homogenizer LD-LMHH-A10

## Overview

Mini Hand-held Homogenizer LX140MHH is a hand-held homogenizer with an easy to dismantle grinding pestle. Equipped with adjustable speed for user control.

## Features:

- Adjustable speed
- Efficient and effective micro-DC motor for grinding trace tissues
- Reuseable grinding pestle with convenient dismantling

## Applications:

Used for grinding various plant tissue sample such as root, stem, leaf, flower, fruit, seed, etc and for grinding various animal tissue sample such as cerebrum, heart, lung, stomach, liver, thymus, kidney, intestine, lymph nodules, muscle, skeleton, etc. Suitable for homogenate.

## Specifications:

<b>Weight</b>	0.18 kg
<b>Battery</b>	AA battery x 2
<b>Dimensions</b>	143 x 25 x 143 mm
<b>Speed range</b>	18000 rpm
<b>Working noise</b>	55 dB
<b>Battery charging time</b>	> 2 hours
<b>Battery non-loaded working time</b>	2 hours