



Cube Ice Makers LD-LCIM-B10

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Cube Ice Makers LD-LCIM-B10

Overview

Cube Ice maker series comes with insulated ice storage bin and door which minimizes air condensation and water formation on equipment surface. It produces steady and similar shape ice cubes with high transparency. It serves safe, hygienic, odorless and best tasting cube ices which can be used for all types of functions. The machine is configured using high quality material and advanced technology at par with ISO and CE standards. It has integrated design for maximum energy saving and water conservation.

Features:

- Corrosion resistant stainless steel is used
- R404A refrigerant is used as refrigerant
- User friendly technology
- Overload protection
- Constructed with stainless steel and rigid plastic cabinets and other high-quality components
- Highly efficient, economical operation
- Fully automatic control, indicator system to indicate of ice and water
- Produces edible ice using potable water
- Ambient temp. 20 Â°C, Inlet water temp. 10 Â°C

Applications:

Cube ice makers are widely used in restaurants and other food service operations for cooling water and beverages. They are ideal for bars to enhance the appearance of alcoholic drinks due to clarity of ice. Often used in laboratories, hotels, marine fisheries, chemical industries, pharmaceutical industries etc.

Specifications:

Weight	230 kgs
Capacity	680 kgs / 24 hrs
Dry Filter	Danfoss(Denmark)
Cooling Fan	EBM(Germany)
Cooling Mode	Air/Water Cooling
Certification	CE
Exterior Finish	304/2B Stainless steel
Compressor Brand	Tecumseh (France)
Standard Voltage	1P / 220 V / 50 Hz
Power Consumption	2800 W
Refrigeration System	R404A

Storage Bin Capacity	460 kgs
Dimension (W x D x H)	840 x 1290 x 830 mm; 1300 x 950 x 1230 mm (2.39CBM)