



## **Multi-tube Vortex Mixer LD-LMVM-A12**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Multi-tube Vortex Mixer LD-LMVM-A12

## Overview

Multi-tube Vortex Mixer LD-LMVM-A12 has a speed range of 500 rpm to 3000 rpm. It uses a brushless DC motor for reliable and smooth operation. The mixer offers touch mode and continuous mode for operation and features an LED display panel for precise control of speed. Its intelligent microprocessor control ensures stable and accurate operation. Our mixer can operate for a maximum of 50 samples at once.

## Features:

- Brushless DC motor
- Uniform acceleration
- Stepless speed regulation
- Wide speed range
- Efficient mixing

## Applications:

Our Multi-tube Vortex Mixer LD-LMVM-A12 is widely used in biological research, clinical diagnostics, pharmaceuticals, and industrial applications for the efficient mixing of multiple samples.

## Specifications:

<b>Racks</b>	Interchangeable racks
<b>Max Load</b>	4.5 kg
<b>Amplitude</b>	4 mm
<b>Dimension</b>	420 × 280 × 335 mm
<b>Speed Range</b>	500 rpm to 3000 rpm
<b>Time Setting</b>	1 min to 24 h
<b>Mixing Method</b>	Horizontal swing
<b>Operation Mode</b>	Touch and continuous mode
<b>Power Consumption</b>	100 W