



## **Platelet Shaker LD-LPS-A13**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Platelet Shaker LD-LPS-A13

## Overview

Platelet shakers are made of high-grade steel with stainless steel shelves for easy cleaning and durability. Supplied with multiple shelves to accommodate 42-84 number of platelet bags. The removable open mesh design provides good air circulation. They are reliable, provides high quality performance and are complied using ISO 9001 & ISO 13485. They are designed to meet the needs of blood banks and clinical laboratories.

## Features:

- Open mesh design provides good air circulation
- Stainless steel shelves can be easily cleaned
- Stationary base drives has storage frame size of 38mm
- Ball bearing drive motor and internal fan for durable and reliable operations
- Robust, practical, and durable design for your economic choice

## Applications:

Platelet shakers are widely used in blood banks and clinical research laboratories.

## Specifications:

<b>Capacity</b>	48 bags
<b>Power Supply</b>	220 V, 60 Hz
<b>Overall Dimension</b>	790 x 255 x 290 mm
<b>Shaking Speed Range</b>	0 - 9 rpm