



**415 L Flammable Storage Cabinet LD-LFSC-  
A15**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# 415 L Flammable Storage Cabinet LD-LFSC-A15

## Overview

415 L Flammable Storage Cabinet LFSC-A15 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Features:

- Capacity - 415 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelf
- Adjustable spill retaining shelves

## Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specifications:

<b>Weight</b>	293 kgs
<b>Capacity</b>	415 L
<b>Door type</b>	Dual door manual
<b>Adjustable shelf</b>	2
<b>Dimension ( H x W x D )</b>	1650 x 1500 x 860 mm
<b>Shelves loading capacity</b>	100 kgs