



**38 L Flammable Storage Cabinet LD-LFSC-  
B11**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# 38 L Flammable Storage Cabinet LD-LFSC-B11

## Overview

38 L Flammable Storage Cabinet LFSC-B11 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Features:

- Capacity - 38 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

## Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specifications:

<b>Weight</b>	56 kgs
<b>Capacity</b>	38 L
<b>Door type</b>	Single door manual
<b>Adjustable shelf</b>	1
<b>Dimension ( H x W x D )</b>	640 x 590 x 600 mm
<b>Shelves loading capacity</b>	50 kgs