



**170 L Flammable Storage Cabinet LD-LFSC-  
C10**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# 170 L Flammable Storage Cabinet LD-LFSC-C10

## Overview

170 L Flammable Storage Cabinet LFSC-C10 is a dual walled self-closing fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Features:

- Capacity - 170 L
- Dual walled self-closing
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

## Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specifications:

|                                 |                        |
|---------------------------------|------------------------|
| <b>Weight</b>                   | 148 kgs                |
| <b>Capacity</b>                 | 170 L                  |
| <b>Door type</b>                | Dual door self-closing |
| <b>Adjustable shelf</b>         | 2                      |
| <b>Dimension ( H x W x D )</b>  | 1650 x 1090 x 460 mm   |
| <b>Shelves loading capacity</b> | 100 kgs                |