



170 L Toxic Cabinet LD-LTXC-A10

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

170 L Toxic Cabinet LD-LTXC-A10

Overview

170 L Toxic Cabinet LTXC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 170 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

Applications:

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

Specifications:

| | |
|---------------------------------|----------------------|
| Weight | 148 kgs |
| Capacity | 170 L |
| Door type | Dual door manual |
| Adjustable shelf | 2 |
| Dimension (H x W x D) | 1650 x 1090 x 460 mm |
| Shelves loading capacity | 100 kgs |