



**114 L Weak Acid and Alkali Cabinet LD-  
LWAC-A10**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# 114 L Weak Acid and Alkali Cabinet LD-LWAC-A10

## Overview

114 L Weak Acid and Alkali Cabinet LWAC-A10 is a dual walled acid and alkali resistance cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

## Features:

- Capacity - 114 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- One adjustable shelf

## Applications:

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

## Specifications:

<b>Weight</b>	117 kgs
<b>Capacity</b>	114 L
<b>Door type</b>	Dual door manual
<b>Adjustable shelf</b>	1
<b>Dimension ( W x D x H )</b>	1090 x 460 x 1120 mm