



**227 L Weak Acid and Alkali Cabinet LD-
LWAC-A12**

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

227 L Weak Acid and Alkali Cabinet LD-LWAC-A12

Overview

227 L Weak Acid and Alkali Cabinet LWAC-A12 is a dual walled manual acid and alkali cabinet with 50 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers. Inner and outer epoxy resin coating provides durability and chemical resistance preventing corrosion and humidity.

Features:

- Capacity - 227 L
- Dual walled manual
- 50 mm insulating air space
- Three point latch and manual lock
- Two adjustable shelves

Applications:

Used to store all kinds of corrosive material or weak acid and alkali, for example: acetic acid, Carbonic acid, phosphoric acid, silicic acid, etc.

Specifications:

| | |
|--------------------------------|---------------------|
| Weight | 173 kgs |
| Capacity | 227 L |
| Door type | Dual door manual |
| Adjustable shelf | 2 |
| Dimension (W x D x H) | 1650 x 860 x 860 mm |