



## **Slide Storage Cabinet LD-LSLC-A10**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Slide Storage Cabinet LD-LSLC-A10

## Overview

Slide Storage Cabinet LD-LSLC-A10 is metallic standalone unit, comes with 18 drawers. Slide cabinet is used to store glass slides, Drawers is consisting of ABS slide plate and cabinet material is cold rolled steel plate and surface is treated with Degree oil, Table regulator, Zinc Parkerizing, Passivation, Powder spraying which is resistant to most chemical aggression.

## Features:

- 4563 or 9072 glass slides can be stored
- 18 Drawers Cabinet
- Slot is made by special mold
- Corrosion resistant material
- Large storage capacity

## Applications:

Slide storage cabinet is used for efficient and organized storage of glass slides across histology, pathology, forensics, research centers etc.

## Specifications:

<b>Capacity</b>	18 Drawers (3 layers x 6 grid, plus a base) Loading: 4536 or 9072 pcs glass slides
<b>Mark groove</b>	Integration forming
<b>Gross weight</b>	27 kgs, 3pcs; 31 kgs, 1pc
<b>Slot material</b>	Made by special mold
<b>Drawer material</b>	Specialized slide drawer, ABS slide plate
<b>Cabinet material</b>	Cold-rolled steel plate
<b>Package Size (W x D x H)</b>	530 x 570 x 440 mm, 3pcs; 530 x 570 x 565mm, 1pc
<b>External Size (W x D x H)</b>	450 x 478 x 1295 mm