



**UV Pass Box (without HEPA filters) LD-LUP-B11**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# UV Pass Box (without HEPA filters) LD-LUP-B11

## Overview

UV Pass Box (without HEPA filters) LUP-B11 is a static type UV pass box equipped UV germicidal lamp, both side electrical interlocked doors, and high quality see through glass window. Pass box is made up of durable high quality 304 grade stainless steel, and updated quality parts to meet the demands of controlled clean room environment operations.

## Features:

- Indicator lights to know the door lock/unlock condition
- Electronically interlocked door type
- Red press button for on-off operations
- Certified with CE and ISO standards

## Applications:

Used to transfer articles, materials between different levels of clean room across pharmaceutical, food, life science and other research, quality inspection industries.

## Specifications:

<b>UV Lamp</b>	30 W × 1
<b>Material</b>	304 Stainless steel
<b>Gross Weight</b>	40 kg
<b>Power Supply</b>	220 V, 50 Hz
<b>Door Interlock Type</b>	Electronical
<b>External Dimension (W × D × H)</b>	760 × 670 × 680 mm
<b>Internal Dimension (W × D × H)</b>	600 × 600 × 600 mm
<b>Packaging Dimension (W × D × H)</b>	800 × 800 × 880 mm