



UV Pass Box (without HEPA filters) LD-LUP-B12

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

UV Pass Box (without HEPA filters) LD-LUP-B12

Overview

UV Pass Box (without HEPA filters) LUP-B12 is a static type UV pass box equipped UV germicidal lamp, both side electrical interlocked doors, and high quality see through glass window. Pass box is made up of durable high quality 304 grade stainless steel, and updated quality parts to meet the demands of controlled clean room environment operations.

Features:

- Indicator lights to know the door lock/unlock condition
- Electronically interlocked door type
- Red press button for on-off operations
- Certified with CE and ISO standards

Applications:

Used to transfer articles, materials between different levels of clean room across pharmaceutical, food, life science and other research, quality inspection industries.

Specifications:

UV Lamp	30 W × 1
Material	304 Stainless steel
Gross Weight	50 kg
Power Supply	220 V, 50 Hz
Door Interlock Type	Electronical
External Dimension (W × D × H)	860 × 770 × 780 mm
Internal Dimension (W × D × H)	700 × 700 × 700 mm
Packaging Dimension (W × D × H)	900 × 900 × 980 mm