



UV Pass Box (without HEPA filters) LD-LUP-B15

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

UV Pass Box (without HEPA filters) LD-LUP-B15

Overview

UV Pass Box (without HEPA filters) LUP-B15 is a static type UV pass box equipped UV germicidal lamp, both side electrical interlocked doors, and high quality see through glass window. Pass box is made up of durable high quality 304 grade stainless steel, and updated quality parts to meet the demands of controlled clean room environment operations.

Features:

- Indicator lights to know the door lock/unlock condition
- Electronically interlocked door type
- Red press button for on-off operations
- Certified with CE and ISO standards

Applications:

Used to transfer articles, materials between different levels of clean room across pharmaceutical, food, life science and other research, quality inspection industries.

Specifications:

UV Lamp	8 to 11 W × 1
Material	304 Stainless steel
Net Weight	55 kg
Power Supply	220 V, 50 Hz
Door Interlock Type	Electronical
External Dimension (W × D × H)	1140 × 1060 × 1080 mm
Internal Dimension (W × D × H)	1000 × 1000 × 1000 mm
Packaging Dimension (W × D × H)	1190 × 1110 × 1180 mm