



## **Mechanical Heating Ultrasonic Cleaner LD-LMUC-A17**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Mechanical Heating Ultrasonic Cleaner LD-LMUC-A17

## Overview

Mechanical Heating Ultrasonic Cleaner LMUC-A17 delivers high efficiency functionality and features. This compact model with 22 L tank capacity and 480 W ultrasonic power provides deep cleansing action removing contamination even from blind holes and internal surfaces. Its key feature is the variance of using it at 28 KHz for gross rinse to 40 KHz for intensive rinse Equipped with degassing and gettering function for removal of dissolved gasses and efficient cleansing of oil and wax

## Features:

- 22 L capacity tank
- Integrated mechanical circuit for timer and heater
- Temperature range : 20 °C ~ 80 °C ( mechanical )
- Cleaning duration : 0 ~ 20 mins
- Efficient cleaning can be done by distilled water
- Optimum durability under prolonged and heavy usage
- Equipped with degassing and gettering function

## Applications:

Ideal for cleaning irregular and complex instruments, lenses and other optical parts, dental and surgical instruments, electronics. It is also used in extraction and to degas certain analyzed stuff in laboratories

## Specifications:

|                                     |                    |
|-------------------------------------|--------------------|
| <b>Voltage</b>                      | AC 220 V / 50 Hz   |
| <b>Capacity</b>                     | 22 L               |
| <b>Control mode</b>                 | Mechanical         |
| <b>Heater power</b>                 | 500 W              |
| <b>Temperature set</b>              | 20 °C ~ 80 °C      |
| <b>Ultrasonic power</b>             | 480 W              |
| <b>Construction material</b>        | Stainless Steel    |
| <b>Tank size ( L x W x H )</b>      | 500 x 300 x 150 mm |
| <b>Machine size ( L x W x H )</b>   | 500 x 330 x 330 mm |
| <b>Digital timer ( mechanical )</b> | 0 ~ 20 mins        |