



## **Mini Vacuum Pump LD-LMVP-A10**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Mini Vacuum Pump LD-LMVP-A10

## Overview

Mini vacuum pump LMVP-A10 uses big preface torque configuration and is particularly designed for low temperature working environment. It portable and requires less storage place. This mini vacuum pump has space-saving design and is ideal for systems integration.

## Features:

- Compact, reliable and powerful
- Especially designed for low temperature working environment
- High pumping speeds with small dimensions
- Eco-friendly exhaust gas free of oil mists
- Quiet, low-vibration operation
- Easy to maintain

## Applications:

It can be used for vacuum packaging, vacuum freezing, gas analysis and other vacuum evacuations.

## Specifications:

<b>Power</b>	90 / 120 W
<b>Weight</b>	6.5kg
<b>Dimension</b>	243 x 114 x 207 mm
<b>Oil Capacity</b>	220 ml
<b>Power Supply</b>	110 V, 60 Hz / 220 V, 50 Hz
<b>Pumping Speed</b>	28.31 to 22.65 L / min
<b>Rotating Speed</b>	1720 / 1440 rpm
<b>Ultimate Pressure</b>	5 Pa