



## **Adherent Cell Culture Flask LD-ACF102L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Adherent Cell Culture Flask LD-ACF102L

## Features:

- Material: High Transparency Polystyrene
- Tissue culture treated
- Sterilization Level is  $10^{-6}$
- Available in both vent cap and plug seal cap
- Gamma ray sterilized
- DNase, RNase and pyrogenic free

## Specifications:

<b>Cap</b>	Plug
<b>Capacity</b>	250 mL
<b>Culture area</b>	$75 \text{ cm}^2$
<b>Quantity per case</b>	200