



## **Suspension Culture Flask LD-SCF103L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Suspension Culture Flask LD-SCF103L

## Features:

- Material: High Transparency Polystyrene
- Sterilization Level is 10<sup>-6</sup>
- Available in both vent cap and plug seal cap
- Gamma ray sterilized
- DNase, RNase and pyrogenic free

## Specifications:

<b>Cap</b>	Vent
<b>Capacity</b>	250 mL
<b>Culture area</b>	75 cm <sup>2</sup>
<b>Quantity per case</b>	200