



Suspension Culture Plate LD-SCP105L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Suspension Culture Plate LD-SCP105L

Features:

- Material: Polystyrene
- Sterilization level is 10^{-6}
- Number codes for easy identification
- Flat surface bottom
- Maximized growth area
- Raised rims on wells to avoid cross contamination
- DNase, RNase and pyrogenic free

Specifications:

Packaging	50 pcs/bag
No of Well	24
Culture area	1.91 cm^2
Packaging Type	Bag
Quantity per case	200