



## **Cell Scraper LD-CSC100L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Cell Scraper LD-CSC100L

## Features:

- Material: Medical grade Polypropylene
- Blade Material: Thermoplastic elastomers
- Gamma ray sterilized
- DNase, RNase and pyrogenic free

## Specifications:

<b>Packaging</b>	1 pcs/bag
<b>Blade Length</b>	1.8 cm
<b>Handle Length</b>	25 cm
<b>Quantity per case</b>	50