



Polyoxymethylene Embedding Cassette LD-PEC109L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polyoxymethylene Embedding Cassette LD-PEC109L

Features:

- Material: Polyoxymethylene
- Cassette writing face is at a 30° and 45° angle
- Cover and base can be disassembled and reassembled with ease
- Resistant to solvents

Specifications:

| | |
|--------------------------|---------------------------|
| Type | 'O" Type Ring |
| Color | Pink |
| Lid Type | No Lid |
| Packaging | 250 pcs/Bag, 20 Bags/Case |
| Quantity Per Case | 5000 |