



## **ABS Centrifuge Tube Rack LD-CTR102L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# ABS Centrifuge Tube Rack LD-CTR102L

## Features:

- Temperature variance range (-80 °C to 121 °C)
- Material: Acrylonitrile Butadiene Styrene
- Autoclavable
- Availability of color customization

## Specifications:

<b>Packaging</b>	Plastic Bag
<b>No. of wells</b>	20 + 30
<b>Tube Capacity</b>	15 mL & 50 mL (Multi use )
<b>Quantity per case</b>	50