



Microcentrifuge Tubes LD-MCT101L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Microcentrifuge Tubes LD-MCT101L

Features:

- Maximum RCF 28000 x g
- Working Temperature is -90 Â°C to 121 Â°C
- Conical bottom design
- Autoclavable temperature is +121 Â°C
- Material: Polypropylene
- Frosted marking area on tube and cap
- Attached cap

Specifications:

Type	Ultra Low Retention
Capacity	0.5/0.65 mL
Packaging	1000 pcs/Bag
Length of tube	30 mm
Quantity per Case	10000
Diameter of the tube	7.2 mm