



## **ABS Cryovials Holder LD-CVH100L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# ABS Cryovials Holder LD-CVH100L

## Features:

- Specially designed for cryovials
- Autoclavable
- Chemically resistant to alcohols and mild organic solvents
- Temperature range is stable from -196 °C to 121 °C

## Specifications:

<b>Capacity</b>	2.0 mL
<b>Diameter</b>	13 mm
<b>Material</b>	Acrylonitrile butadiene styrene
<b>Packaging</b>	10 pcs/Bag
<b>Number of Wells</b>	25
<b>Quantity per Case</b>	100
<b>Dimension ( L*W*H )</b>	76*76*54 mm